

## Anti-Phospho-ATXN3-Ser256 antibody (STJA0009980)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                     |                                                                                                           |
|---------------------|-----------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                |
| Concentration       | 1 mg/mL                                                                                                   |
| Conjugation         | Unconjugated                                                                                              |
| Purification        | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Dilution Range      | WB 1:500-2000                                                                                             |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                   |
| Isotype             | IgG                                                                                                       |
| Molecular Weight    | Calculated: 40kD                                                                                          |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.            |

### TARGET INFORMATION

|             |                                                                                           |
|-------------|-------------------------------------------------------------------------------------------|
| Gene ID     | <a href="#">4287</a>                                                                      |
| Gene Symbol | <a href="#">ATXN3</a>                                                                     |
| UniProt ID  | <a href="#">ATX3_HUMAN</a>                                                                |
| Specificity | This antibody detects endogenous levels of Ataxin-3 (Phospho-Ser256) at Human, Mouse, Rat |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Western Blot analysis of various, using primary antibody at 1:1000 dilution. Secondary antibody (STJS000791) was diluted at 1:10000