

## Anti-Phospho-RBL2-Thr642 antibody (STJA0010062)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC/WB             |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                                           |
|---------------------|-----------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                |
| Concentration       | 1 mg/mL                                                                                                   |
| Conjugation         | Unconjugated                                                                                              |
| Purification        | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Dilution Range      | WB 1:500-2000<br>IHC 1:50-200                                                                             |
| Formulation         | Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.                                             |
| Isotype             | IgG                                                                                                       |
| Molecular Weight    | Calculated: 125kD                                                                                         |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.            |

### TARGET INFORMATION

|             |                                                                                                  |
|-------------|--------------------------------------------------------------------------------------------------|
| Gene ID     | <a href="#">5934</a>                                                                             |
| Gene Symbol | <a href="#">RBL2</a>                                                                             |
| UniProt ID  | <a href="#">RBL2_HUMAN</a>                                                                       |
| Specificity | This antibody detects endogenous levels of RBL2 (Phospho Thr642) Rabbit pAb at Human, Mouse, Rat |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western Blot analysis of 1 HeLa cell, 22 Serum-free treated, using primary antibody at 1:1000 dilution. Secondary antibody (STJS000791) was diluted at 1:10000