

Anti-V5-Tag antibody [SM0385] (STJA0010385)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/ICC/IF
Host / Source	Mouse
Reactivity	V5-Tag

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1k
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Immunofluorescent analysis of 4% paraformaldehyde-fixed transfected HEK-293 cells using V5 Tag at dilution of 1:1000



Immunofluorescent analysis of 4% paraformaldehyde-fixed transfected HEK-293 cells using V5 Tag at dilution of 1:1000