

Anti-CD19 antibody [SM0393] (STJA0010393)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H) 1:500-1:1000 IHC/IF (H) 1:500-1:1500
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG2ak
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	930
Gene Symbol	CD19
UniProt ID	CD19_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot analysis of CD19 at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded human tonsil using CD19 at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded human lymphoma using CD19 at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded human lung squamous cell carcinoma using CD19 at dilution of 1:1000