

Anti-Claudin18.2 antibody [SM0395] (STJA0010395)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H, R) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG2ak
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Immunohistochemistry analysis of paraffin embedded rat stomach using claudin18.2 at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded human stomach using claudin18.2 at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded human gastric carcinoma using claudin18.2 at dilution of 1:500