

Anti-GFAP antibody [SM0405] (STJA0010405)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (M, R) 1:500-1:5000 IHC/IF (H, M, R) 1:500-1:2000/1:500-1:2000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	14580
Gene Symbol	Gfap
UniProt ID	GFAP_MOUSE

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




Immunohistochemistry analysis of GFAP. Sample: mouse brain, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, 98C, 20 min. Blocking buffer: 3% BSA in PBS, RT, 30min. Primary antibody: 1:600, 4C overnight. Secondary antibody: HRP Goat Anti-mouse IgG, 1:200 RT 1h.




Immunohistochemistry analysis of GFAP. Sample: Rat brain, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, 98C, 20 min. Blocking buffer: 3% BSA in PBS, RT, 30min. Primary antibody: 1:600, 4C overnight. Secondary antibody: HRP Goat Anti-mouse IgG, 1:200 RT 1h.



IF analysis of GFAP. Sample: mouse brain, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, 98C, 20 min. Blocking buffer: 3% BSA in PBS, RT, 30min. Primary antibody: 1:1000, 4C overnight. Secondary antibody: Cy3 conjugated Goat Anti-mouse IgG, 1:200 RT 1h.



IF analysis of GFAP. Sample: Rat brain, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, 98C, 20 min. Blocking buffer: 3% BSA in PBS, RT, 30min. Primary antibody: 1:1000, 4C overnight. Secondary antibody: Cy3 conjugated Goat Anti-mouse IgG, 1:200 RT 1h.



IF analysis of GFAP. Sample: Rat brain, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, 98C, 20 min. Blocking buffer: 3% BSA in PBS, RT, 30min. Primary antibody: 1:1000, 4C overnight. Secondary antibody: Cy3 conjugated Goat Anti-mouse IgG, 1:200 RT 1h.