

Anti-GFAP antibody [SM0407] (STJA0010407)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Mouse/Rat

PRODUCT PROPERTIES



Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (M, R) 1:500-1:1000 IHC/IF (H, M, R, Mk) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG2bk
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2670
Gene Symbol	GFAP
UniProt ID	GFAP_HUMAN

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




Immunohistochemistry analysis of paraffin embedded rat brain using GFAP at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded mouse brain using GFAP at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded human glioma using GFAP at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded human oligodendroglioma using GFAP at dilution of 1:1000

