

Anti-CD10 antibody [SM0419] (STJA0010419)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (M) 1:1000-1:3000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG2bk
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4311
Gene Symbol	MME
UniProt ID	NEP_HUMAN

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Immunohistochemistry analysis of paraffin embedded mouse prostate using CD10 at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded mouse small intestine using CD10 at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded mouse kidney using CD10 at dilution of 1:1000