

## Anti-Ki67 antibody [SM0421] (STJA0010421)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/IF          |
| Host / Source | Mouse              |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES




|                     |                                                                          |
|---------------------|--------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                               |
| Conjugation         | Unconjugated                                                             |
| Purification        | Affinity purification                                                    |
| Dilution Range      | WB (H, M, R) 1:500-1:1000<br>IHC/IF (H) 1:500-1:1500                     |
| Formulation         | PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.             |
| Isotype             | IgG1k                                                                    |
| Storage Instruction | Store at -20C for up to one year, and avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">4288</a>       |
| Gene Symbol | <a href="#">MKI67</a>      |
| UniProt ID  | <a href="#">KI67_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin embedded human colon using ki67 at dilution of 1:1000

Immunohistochemistry analysis of paraffin embedded human tonsil using ki67 at dilution of 1:1000