

Anti-RFP antibody [SM0494] (STJA0010494)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	ICC/IF
Host / Source	Mouse
Reactivity	RFP

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ICC/IF (Species independent) 1:500-1:1500
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG2ak
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Immunofluorescent analysis of 4% paraformaldehyde-fixed transfected HEK-293 cells using RFP Tag at dilution of 1:500