

Anti-Cytokeratin 6 antibody [SM0507] (STJA0010507)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H, M, R) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1k
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3853
Gene Symbol	KRT6A
UniProt ID	K2C6A_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Immunohistochemistry analysis of paraffin embedded rat skin using CK6 at dilution of 1:500




Immunohistochemistry analysis of paraffin embedded human esophagus using CK6 at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded mouse skin using CK6 at dilution of 1:500



Immunohistochemistry analysis of paraffin embedded human cervical carcinoma using CK6 at dilution of 1:500



Immunohistochemistry analysis of paraffin embedded mouse epithelial of oral cavity using CK6 at dilution of 1:500