

Anti-Recombinant-Progesterone receptor antibody [SR0519] (STJA0010519)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	WB (H) 1:500-1:1000 IHC/IF (H, M, R) 1:200-1:400
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5241
Gene Symbol	PGR
UniProt ID	PRGR_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

