

Anti-Recombinant-Androgen receptor antibody [SR0520] (STJA0010520)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	IHC/IF (H, M, R) 1:200-1:400
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	367
Gene Symbol	AR
UniProt ID	ANDR_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---