

**Anti-Recombinant-SOX9 antibody [SR0521] (STJA0010521)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC/IF             |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

**PRODUCT PROPERTIES**

|                     |                                                                         |
|---------------------|-------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                              |
| Conjugation         | Unconjugated                                                            |
| Purification        | Recombinant Expression and Affinity purified                            |
| Dilution Range      | IHC/IF (H, M, R) 1:200-1:400                                            |
| Formulation         | PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA                          |
| Storage Instruction | Store at-20C for up to one year, and avoid repeated freeze-thaw cycles. |

**TARGET INFORMATION**

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">6662</a>       |
| Gene Symbol | <a href="#">SOX9</a>       |
| UniProt ID  | <a href="#">SOX9_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|