

Anti-Recombinant-Phospho-Granzyme B antibody [SR0537] (STJA0010537)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (R) 1:250-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgGk
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3002
Gene Symbol	GZMB
UniProt ID	GRAB_HUMAN

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




Immunohistochemistry analysis of Granzyme B. Sample: Rat thymus, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, pressure cooker 2 min. Primary antibody: 1:800, 4C overnight. Secondary antibody: S-vision poly-HRP conjugated Goat Anti-Rabbit IgG, Ready to use, RT, 20min.




Immunohistochemistry analysis of Granzyme B. Sample: Rat thymus, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, pressure cooker 2 min. Primary antibody: 1:800, 4C overnight. Secondary antibody: S-vision poly-HRP conjugated Goat Anti-Rabbit IgG, Ready to use, RT, 20min.




Immunohistochemistry analysis of Granzyme B. Sample: Rat thymus, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, pressure cooker 2 min. Primary antibody: 1:800, 4C overnight. Secondary antibody: S-vision poly-HRP conjugated Goat Anti-Rabbit IgG, Ready to use, RT, 20min.



Immunohistochemistry analysis of Granzyme B. Sample: Rat thymus, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, pressure cooker 2 min. Primary antibody: 1:800, 4C overnight. Secondary antibody: S-vision poly-HRP conjugated Goat Anti-Rabbit IgG, Ready to use, RT, 20min.



Immunohistochemistry analysis of Granzyme B. Sample: Rat thymus, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, pressure cooker 2 min. Primary antibody: 1:1000, 4C overnight. Secondary antibody: S-vision poly-HRP conjugated Goat Anti-Rabbit IgG, Ready to use, RT, 20min.



Immunohistochemistry analysis of Granzyme B. Sample: Rat thymus, 4% PFA 12-24h. Antigen retrieval: Citrate buffer, pressure cooker 2 min. Primary antibody: 1:1000, 4C overnight. Secondary antibody: S-vision poly-HRP conjugated Goat Anti-Rabbit IgG, Ready to use, RT, 20min.