

## Anti-Recombinant-Phospho-JNK-T183/Y185 antibody [SR0546] (STJA0010546)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                          |
|---------------------|--------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                               |
| Conjugation         | Unconjugated                                                             |
| Purification        | Affinity purification                                                    |
| Dilution Range      | WB (H, M, R) 1:300-1:600                                                 |
| Formulation         | PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.             |
| Isotype             | IgG                                                                      |
| Storage Instruction | Store at -20C for up to one year, and avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|         |                                                                                                                                                                                                                   |
|---------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID | <a href="#">5601</a><br><a href="#">5602</a><br><a href="#">5599</a><br><a href="#">MAPK9</a><br><a href="#">MAPK10</a><br><a href="#">MK09_HUMAN</a><br><a href="#">MK10_HUMAN</a><br><a href="#">MK08_HUMAN</a> |
|---------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Western blot analysis of P-JNK at dilution of 1:600