

Anti-Recombinant-GRP78 BiP antibody [SM0567] (STJA0010567)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC/IF
	Host / Source	Mouse
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB (H, M, R) 1:500-1:1000 IHC/IF (H, M, R) 1:500-1:2000
	Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
	Isotype	IgG1
	Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	14828
	Gene Symbol	Hspa5
	UniProt ID	BIP_MOUSE
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.