

## Anti-Recombinant-Cytokeratin 1 antibody [SM0571] (STJA0010571)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Purification        | Affinity purification   |
| Dilution Range      | WB (H) 1:500-1:1000   |
| Formulation         | PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.            |
| Isotype             | IgG1  |
| Storage Instruction | Store at-20C for up to one year, and avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">3848</a>       |
| Gene Symbol | <a href="#">KRT1</a>       |
| UniProt ID  | <a href="#">K2C1_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|