

**Anti-Recombinant-Hexokinase II antibody [SM0575] (STJA0010575)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Mouse              |
| Reactivity    | Human/Mouse/Rat    |

**PRODUCT PROPERTIES**

|                     |                                                                          |
|---------------------|--------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                               |
| Conjugation         | Unconjugated                                                             |
| Purification        | Affinity purification                                                    |
| Dilution Range      | WB (H, M, R) 1:500-1:1000                                                |
| Formulation         | PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.             |
| Isotype             | IgG1                                                                     |
| Storage Instruction | Store at -20C for up to one year, and avoid repeated freeze-thaw cycles. |

**TARGET INFORMATION**

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">3099</a>       |
| Gene Symbol | <a href="#">HK2</a>        |
| UniProt ID  | <a href="#">HXK2_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|