

## Anti-Recombinant-Glutamine Synthetase antibody [SM0584] (STJA0010584)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:500-1:1000 IHC/IF (H, M, R) 1:3000-1:6000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2752</a>
Gene Symbol	<a href="#">GLUL</a>
UniProt ID	<a href="#">GLNA_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---