

## Anti-Recombinant-CDX2 antibody [SM0606] (STJA0010606)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/IF          |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                         |
|---------------------|-------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                              |
| Conjugation         | Unconjugated                                                            |
| Purification        | Affinity purification                                                   |
| Dilution Range      | WB (H) 1:500-1:2000<br>IHC/IF (H, M, R) 1:500-1:2000                    |
| Formulation         | PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.            |
| Isotype             | IgG1                                                                    |
| Storage Instruction | Store at-20C for up to one year, and avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">1045</a>       |
| Gene Symbol | <a href="#">CDX2</a>       |
| UniProt ID  | <a href="#">CDX2_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of CDX2. Sample: Protein treated by RIPA Lysis Buffer. Blocking buffer: 3% Nonfat dry milk in TBST, RT, 1h. Primary antibody: 1:2000, 4C overnight. Secondary antibody: HRP conjugated Goat Anti-mouse IgG, 1:3000, RT, 1h.