

## Anti-Recombinant-PGP9.5 antibody [SM0609] (STJA0010609)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">7345</a>
Gene Symbol	<a href="#">UCHL1</a>
UniProt ID	<a href="#">UCHL1_HUMAN</a>

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of PGP9. 5 at dilution of 1:1000