

Anti-Recombinant-Cytokeratin 15 antibody [SM0610] (STJA0010610)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC/IF
	Host / Source	Mouse
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB (H) 1:500-1:1500 IHC/IF (H, R) 1:250-1:1000
	Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
	Isotype	IgG1
	Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	3866
	Gene Symbol	KRT15
	UniProt ID	K1C15_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.