

## Anti-Recombinant-Tyrosine Hydroxylase antibody [SR0614] (STJA0010614)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC/IF             |
| Host / Source | Rabbit             |
| Reactivity    | Mouse/Rat/Monkey   |

### PRODUCT PROPERTIES

|                     |                                                                         |
|---------------------|-------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                              |
| Conjugation         | Unconjugated                                                            |
| Purification        | Affinity purification                                                   |
| Dilution Range      | IHC/IF (M, R, Mk) 1:500-1:1000                                          |
| Formulation         | PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.            |
| Isotype             | IgG                                                                     |
| Storage Instruction | Store at-20C for up to one year, and avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">21823</a>      |
| Gene Symbol | <a href="#">Th</a>         |
| UniProt ID  | <a href="#">TY3H_MOUSE</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|

