

Anti-Recombinant-alpha Tubulin antibody [SR0621] (STJA0010621)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:2000-1:4000 IHC/IF (H, M, R) 1:2000-1:4000 ICC/IF (H, M, R) 1:2000-1:4000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7277
Gene Symbol	TUBA4A
UniProt ID	TBA4A_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of Alpha-Tubulin at dilution of 1:2000



Immunofluorescence analysis of paraffin embedded mouse kidney using Alpha-Tubulin at dilution of 1:1000



Immunocytochemistry analysis of 4% paraformaldehyde-fixed HeLa cells using Alpha-Tubulin at dilution of 1:1000