

Anti-Recombinant-Bcl-2 antibody [SM0676] (STJA0010676)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB
	Host / Source	Mouse
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB (H, M, R) 1:1000-1:2000
	Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
	Isotype	IgG1
	Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	23786
	Gene Symbol	BCL2L13
	UniProt ID	B2L13_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.