

**Anti-Recombinant-HA-Tag antibody [SR0721] {HRP} (STJA0010721)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	HA-Tag

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:1000-1:2000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---