

Anti-Recombinant-MBP-Tag antibody [SM0726] {HRP} (STJA0010726)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	MBP-Tag

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:500-1:1500
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---