

## Anti-Recombinant-COX IV antibody [SM0729] {HRP} (STJA0010729)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB (M, R) 1:1000-1:2000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">12857</a>
Gene Symbol	<a href="#">Cox4i1</a>
UniProt ID	<a href="#">COX41_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---