

Anti-APP antibody (66-82 aa) [R09-3F7] (STJA0010881)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	351
Gene Symbol	APP
UniProt ID	A4_HUMAN
Immunogen Region	66-82 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Amyloid Precursor Protein antibody western blot mouse brain STJA0010881



Amyloid Precursor Protein antibody western blot rat brain STJA0010881



Amyloid Precursor Protein antibody immunohistochemistry human breast cancer STJA0010881



Amyloid Precursor Protein antibody immunocytochemistry HeLa cells STJA0010881