

Anti-TUBA4A antibody (400-448 aa) [R02-2H3] (STJA0010943)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Hamster

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	7277
Gene Symbol	TUBA4A
UniProt ID	TBA4A_HUMAN
Immunogen Region	400-448 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




alpha tubulin antibody western blot mouse brain 3T3 STJA0010943



alpha tubulin antibody western blot HeLa CHO-K1 rat brain STJA0010943



alpha Tubulin antibody immunohistochemistry human tonsil STJA0010943



alpha tubulin antibody immunocytochemistry HeLa cells STJA0010943