

Anti-ATP5F1B antibody (151-164 aa) [R01-4G1] (STJA0011009)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Monkey

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	506
Gene Symbol	ATP5F1B
UniProt ID	ATPB_HUMAN
Immunogen Region	151-164 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ATPB antibody western blot mouse rat heart lysates STJA0011009



ATPB antibody western blot Jurkat HeLa COS7 Raji STJA0011009



ATPB antibody immunohistochemistry human tonsil STJA0011009



ATPB antibody immunocytochemistry HeLa cells STJA0011009