

Anti-TUBB4A antibody (430-444 aa) [R06-6A5] (STJA0011034)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	10382
Gene Symbol	TUBB4A
UniProt ID	TBB4A_HUMAN
Immunogen Region	430-444 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




beta IV Tubulin antibody western blot mouse brain STJA0011034



beta IV Tubulin antibody western blot K562 rat brain C6 HeLa STJA0011034



Beta-IV tubulin antibody immunocytochemistry HeLa cells STJA0011034



beta IV Tubulin antibody immunohistochemistry human brain STJA0011034