

## Anti-Phospho-CREB1-Ser133 antibody (26-140 aa) [R07-3C4] (STJA0011169)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P/IP        |
| Host / Source | Rabbit             |
| Reactivity    | Human/Rat          |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.3 mg/mL                                                                                |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100<br>IP 1:20                                             |
| Formulation         | 50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA   |
| Isotype             | IgG                                                                                      |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">1385</a>        |
| Gene Symbol      | <a href="#">CREB1</a>       |
| UniProt ID       | <a href="#">CREB1_HUMAN</a> |
| Immunogen Region | 26-140 aa                   |

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Phospho-CREB Ser133 antibody western blot C6 UV lysates STJA0011169



Phospho-CREB Ser133 antibody western blot CHO-K1 STJA0011169



CREB Phospho-Ser133 antibody immunohistochemistry human lung cancer STJA0011169