

Anti-SRSF1 antibody (16-30 aa) [R07-9G7] (STJA0011952)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	6426
Gene Symbol	SRSF1
UniProt ID	SRSF1_HUMAN
Immunogen Region	16-30 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




SF2 antibody immunocytochemistry 293T cells
STJA0011952



SF2 antibody western blot 3T3 lysates
STJA0011952



SF2 antibody western blot K562 C6 lysates
STJA0011952



SF2 antibody immunohistochemistry human tonsil
STJA0011952