

Anti-POLR2A antibody (1602-1610 aa) [R06-8G7] (STJA0012190)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	5430
Gene Symbol	POLR2A
UniProt ID	RPB1_HUMAN
Immunogen Region	1602-1610 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



RNA Polymerase II Subunit B1 antibody immunocytochemistry CEM STJA0012190



RNA polymerase II CTD repeat YSPTSPS antibody western blot HeLa A549 MCF-7 STJA0012190



RNA polymerase II CTD repeat YSPTSPS antibody immunohistochemistry human breast cancer STJA0012190