

## Anti-EEF1A1 antibody (325-338 aa) [R04-3H4] (STJA0012354)

### GENERAL INFORMATION

|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | WB/IHC-F/IHC-P/ICC/IF |
| Host / Source | Rabbit                |
| Reactivity    | Human/Rat             |

### PRODUCT PROPERTIES

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.3 mg/mL                                                                                |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100<br>IF 1:50-1:200                                       |
| Formulation         | 50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA   |
| Isotype             | IgG                                                                                      |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">1915</a>        |
| Gene Symbol      | <a href="#">EEF1A1</a>      |
| UniProt ID       | <a href="#">EF1A1_HUMAN</a> |
| Immunogen Region | 325-338 aa                  |

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Elongation Factor 1A1 antibody immunocytochemistry HeLa cells STJA0012354



eEF1A1/EFTu antibody western blot U2OS C6 STJA0012354



eEF1A1/EFTu antibody immunohistochemistry human tonsil STJA0012354