

Anti-FXR1 antibody (42-56 aa) [R03-011] (STJA0012382)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	8087
Gene Symbol	FXR1
UniProt ID	FXR1_HUMAN
Immunogen Region	42-56 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




FXR1 antibody immunocytochemistry HeLa cells
STJA0012382



FXR1 antibody western blot PC-12 C2C12
STJA0012382



FXR1 antibody western blot HeLa MCF-7 HepG2
U2OS U251 HL-60 STJA0012382



FXR1 antibody immunohistochemistry human
lung cancer STJA0012382