

Anti-IMPDH2 antibody (140-154 aa) [R06-1B8] (STJA0012430)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	3615
Gene Symbol	IMPDH2
UniProt ID	IMDH2_HUMAN
Immunogen Region	140-154 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




IMPDH2 antibody immunocytochemistry HeLa cells STJA0012430



IMPDH2 antibody western blot 3T3 mouse liver STJA0012430



IMPDH2 antibody western blot HeLa A549 HL-60 U2OS C6 STJA0012430



IMPDH2 antibody immunohistochemistry human breast cancer STJA0012430