

## Anti-NPC2 antibody (20-151 aa) [R07-0C8] (STJA0012496)

### GENERAL INFORMATION

|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | WB/IHC-F/IHC-P/ICC/IF |
| Host / Source | Rabbit                |
| Reactivity    | Human/Mouse/Rat       |

### PRODUCT PROPERTIES


|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 0.3 mg/mL  |
| Conjugation         | Unconjugated   |
| Purification        | Affinity Purified  |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100<br>IF 1:50-1:200                                       |
| Formulation         | 50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA   |
| Isotype             | IgG  |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">10577</a>      |
| Gene Symbol      | <a href="#">NPC2</a>       |
| UniProt ID       | <a href="#">NPC2_HUMAN</a> |
| Immunogen Region | 20-151 aa                  |

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Niemann Pick C2 antibody immunocytochemistry U87-MG STJA0012496



Niemann Pick C2 antibody western blot Raw2647 PC-12 C6 STJA0012496



Niemann Pick C2 antibody western blot HeLa A549 HL-60 STJA0012496



Niemann Pick C2 antibody immunohistochemistry human lung cancer STJA0012496