

Anti-UBE4B antibody (28-40 aa) [R04-6C6] (STJA0012638)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IF 1:50-1:200
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	10277
Gene Symbol	UBE4B
UniProt ID	UBE4B_HUMAN
Immunogen Region	28-40 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



UBE4B antibody western blot PC-12 C6
STJA0012638



UBE4B antibody western blot HeLa A549 HL-60
U2OS U897MG STJA0012638



UBE4B antibody immunocytochemistry MCF-7
cells STJA0012638