

Anti-PRKCB antibody (170-350 aa) [R06-0C2] (STJA0012652)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	5579
Gene Symbol	PRKCB
UniProt ID	KPCB_HUMAN
Immunogen Region	170-350 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




PKC beta 2 antibody immunocytochemistry K562 cells STJA0012652



PKC beta 2 antibody western blot rat brain lysates STJA0012652



PKC beta 2 antibody western blot HeLa A549 U2OS U87-MG STJA0012652



PKC beta 2 antibody immunohistochemistry human tonsil paraffin-embedded STJA0012652