

Anti-RPA1 antibody (317-616 aa) [8C3-D12-H10] (STJA0012829)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF/IP
Host / Source	Mouse
Reactivity	Human/Monkey/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IF 1:50-1:200 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2a
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	6117
Gene Symbol	RPA1
UniProt ID	RFA1_HUMAN
Immunogen Region	317-616 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




RPA70 antibody western blot HeLa A431 MCF-7
COS7 HT-1080 K562 STJA0012829




RPA70 antibody immunocytochemistry HeLa cells
STJA0012829



RPA7 antibody immunoprecipitation HeLa lysates
STJA0012829



RPA70 8C3 antibody western blot HeLa 293T C6
3T3 K562 STJA0012829



RPA70 antibody western blot HUVEC HaCat 3T3
MCF-7 A549 STJA0012829