

## Anti-DDX3X antibody (22-178 aa) [6G8-F4-E3] (STJA0012848)

### GENERAL INFORMATION

|               |                        |
|---------------|------------------------|
| Product Type  | Primary antibodies     |
| Applications  | WB/ICC/IF/IP/ChIP      |
| Host / Source | Mouse                  |
| Reactivity    | Human/Rat/Mouse/Monkey |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 3.41 mg/mL                                                                               |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IF 1:50-1:200<br>IP 1:20<br>ChIP 1:20                                 |
| Formulation         | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.          |
| Isotype             | IgG2a                                                                                    |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">1654</a>        |
| Gene Symbol      | <a href="#">DDX3X</a>       |
| UniProt ID       | <a href="#">DDX3X_HUMAN</a> |
| Immunogen Region | 22-178 aa                   |

### ADDITIONAL INFORMATION


**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



DDX3 antibody western blot HeLa 3T3 C6 COS7 K562 Jurkat STJA0012848



DDX3 antibody immunocytochemistry HeLa STJA0012848



DDX3 antibody immunoprecipitation HeLa lysates STJA0012848