

Anti-RELA antibody (247-296 aa) [3D2-4E9-7A8] (STJA0012897)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1.43 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	5970
Gene Symbol	RELA
UniProt ID	TF65_HUMAN
Immunogen Region	247-296 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




NF-KB p65 antibody western blot T47D 3T3 COS7 C6 HeLa STJA0012897




NF-KB p65 antibody immunofluorescence HeLa cells STJA0012897



NF-KB p65 antibody immunoprecipitation HeLa lysates STJA0012897



NF-KB p65 antibody immunohistochemistry rat brain tissue STJA0012897



NF-KB p65 antibody western blot COS7 C6 T47D STJA0012897