

Anti-SOX2 antibody (1-68 aa) [R74-1A-2] (STJA0016073)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000 IHC 1:100-1:200 IF 1:50-1:200
Formulation	Liquid in PBS, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	6657
Gene Symbol	SOX2
UniProt ID	SOX2_HUMAN
Immunogen Region	1-68 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SOX2 antibody western blot NCCIT F9 lysates
STJA0016073



SOX2 antibody immunohistochemistry human
cervix carcinoma STJA0016073



SOX2 antibody immunofluorescence NCCIT cells
STJA0016073