

Syrian Hamster IgG Isotype Control [SHG-1] {EV500} (STJA0016279)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Syrian Hamster

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	EV500
Purification	Affinity purified
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---