

## Rat IgG1k Isotype Control [HRPN] {EV610} (STJA0016298)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rat                |

### PRODUCT PROPERTIES

|                     |                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                       |
| Conjugation         | EV610                                                                                            |
| Purification        | Affinity purified                                                                                |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG1k                                                                                            |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### ADDITIONAL INFORMATION

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|