

## Anti-CD3 antibody [17A2] {EF700} (STJA0016328)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	EF700
Purification	Affinity purified
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">12500</a> <a href="#">12502</a> <a href="#">12503</a> <a href="#">Cd3d</a> <a href="#">Cd3g</a> <a href="#">CD3D_MOUSE</a> <a href="#">CD3G_MOUSE</a> <a href="#">CD3Z_MOUSE</a>
---------	---

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---