

Anti-CD48 antibody [HM48-1] {EF700} (STJA0016330)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	FC
	Host / Source	Armenian Hamster
	Reactivity	Mouse
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	EF700
	Purification	Affinity purified
	Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
	Isotype	IgG
	Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
TARGET INFORMATION		
	Gene ID	12506
	Gene Symbol	Cd48
	UniProt ID	CD48_MOUSE
ADDITIONAL INFORMATION		
	Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.